

## TestOxidase™ Reagent



- Rapid detection of bacterial cytochrome oxidase
- Stored at room temperature
- Nine months shelf life
- Simple to use
- Up to 400 tests per bottle
- Flexible applications
- Convenient dropper bottle

Code	Product	Size
PL.390	TestOxidase™ Reagent	15 ml

## Spot Indole Reagent



- Rapid detection of tryptophan degradation
- Simple to use
- Clear results
- Flexible applications
- Convenient dropper bottle

Code	Product	Size
PL.391-10	Spot Indole Reagent	10 ml

## Mucolyse Sputum Digestant



- Effective sputum digestant
- Simple and convenient to use
- Economical and safe format
- Improved isolation of potential pathogens
- No effect on morphology, growth of FA standing of potential pathogens
- Two year shelf life
- Simply reconstitutes to 100 ml of active solution

Code	Product	Size
PL.701	Mucolyse Sputum Digestant	10 Vials / Box

Pro-Lab Diagnostics U.K.  
 3 Bassendale Road  
 Bromborough, Wirral, Merseyside  
 CH62 3QL  
 United Kingdom  
 Tel: 0151 353 1613  
 Fax: 0151 353 1614  
 E-mail: [uksupport@pro-lab.com](mailto:uksupport@pro-lab.com)

Pro-Lab Diagnostics Canada  
 20 Mural Street, Unit #4  
 Richmond Hill, ON L4B 1K3  
 Canada  
 Toll Free: 1-800-268-2341  
 Tel: (905) 731-0300  
 Fax: (905) 731-0206  
 E-mail: [support@pro-lab.com](mailto:support@pro-lab.com)

Pro-Lab Diagnostics U.S.A.  
 21 Cypress Blvd., Suite 1070  
 Round Rock, TX, 78665-1034  
 U.S.A.  
 Toll Free: 1-800-522-7740  
 Tel: (512) 832-9145  
 Fax: 1-800-332-0450  
 E-mail: [ussupport@pro-lab.com](mailto:ussupport@pro-lab.com)



[www.pro-lab.com](http://www.pro-lab.com)

## Pro-Loops™ Inoculation Loops and Holders



- Calibrated nichrome wire loops
- Certified nichrome wire loops
- Range of 6 sizes from 1 µl to 10 µl
- High quality nickel plated brass and light weight aluminum holders with insulated handles
- Holders available in 6" and 8" sizes

Code	Product	Volume	Diameter	Size
PL.100	10 µl Nichrome Wire Loops	1/100 ml	5.05 mm	25 Loops
PL.101	5 µl Nichrome Wire Loops	1/200 ml	3.57 mm	25 Loops
PL.102	3.3 µl Nichrome Wire Loops	1/300 ml	2.91 mm	25 Loops
PL.103	2.5 µl Nichrome Wire Loops	1/400 ml	2.52 mm	25 Loops
PL.104	2 µl Nichrome Wire Loops	1/500 ml	2.26 mm	25 Loops
PL.105	1 µl Nichrome Wire Loops	1/1000 ml	1.59 mm	25 Loops
PL.130	10 µl Heavy Gauge Wire Loops	1/100 ml	5.05 mm	25 Loops
PL.131	5 µl Heavy Gauge Wire Loops	1/200 ml	3.57 mm	25 Loops
PL.133	2.5 µl Heavy Gauge Wire Loops	1/400 ml	2.52 mm	25 Loops
PL.106/6-1	Insulated Nickel Plated Brass Holder			6 Inch
PL.106/8-1	Insulated Nickel Plated Brass Holder			8 Inch
PL.107/6-1	Insulated Light Weight Aluminum Holder			6 Inch
PL.107/8-1	Insulated Light Weight Aluminum Holder			8 Inch

## Certified / Calibrated Nichrome Wire Loops

- Manufactured from the highest quality 80/20 nickel chromium wire
- Calibrated using the Evans Blue Dye Method (Cumitech 2A)
- Each pack of 5 loops is presented in a protective case
- Accompanied by the appropriate certificate of calibration

Code	Product	Volume	Diameter	Size
PL.120	10 µl Certified / Calibrated Wire Loops	1/100 ml	5.05 mm	5 Loops
PL.125	1 µl Certified / Calibrated Wire Loops	1/1000 ml	1.59 mm	5 Loops

## Disposable Plastic Inoculation Loops

A range of plastic disposable loops which are calibrated. Each case contains a certificate of calibration.

Code	Product	Colour	Size
PL.W710/20	10 µl Calibrated Plastic Disposable Loops	Red	Case of 3000
PL.W711/20	1 µl Calibrated Plastic Disposable Loops	Blue	Case of 3000

## McFarland Standards



- Hazard-free format
- Two year shelf life
- Stored at room temperature
- Supplied in protective carrying case
- Full range of McFarland Standards available
- Presented in convenient 100mm x 15mm tubes

Code	Product	Size
SD2350	McFarland Standard Set Containing - 1 each of 0.5, 1.0, 2.0, 3.0, 4.0	Each
SD2300	McFarland Standard 0.5	Each
SD2301	McFarland Standard 1.0	Each
SD2302	McFarland Standard 2.0	Each
SD2303	McFarland Standard 3.0	Each
SD2304	McFarland Standard 4.0	Each

## EnviroBead™ Autoclave Deodorants



- Effective autoclave deodorant
- Odours changed and eliminated, not masked
- Simple to use
- Naturally extracted fragrances
- 100 EnviroBead™ capsules per jar
- One EnviroBead™ will deodorize 10 cubic feet of waste

**"Neutralize - Don't mask"**

Code	Product	Size
PL.603/100	Peach EnviroBead™	100 Capsules
PL.604/100	Apple EnviroBead™	100 Capsules
PL.605/100	Lemon EnviroBead™	100 Capsules